

# 'Preserved Flowers' You can make it by yourself !



Rose-Red



Classic-Red



Carbon-Black

● **Mysty** is the processing liquid which operation is simple.  
Please see photographs below, you can easily process it.



### Fresh flower

Please choose your favorite flower.




### Decolorization

Please prepare a carafe with a cover.  
Soak it in solution for decolorization for 24hours.



### Coloration

Soak the flower, which you decolorized in solution of a favorite color for coloration for 24-48 hours.  
You can mix solution of favorite colors.



### Drying

Please dry it at room temperature well.



### Completion

It is completion.  
Please keep it with care not to get damp.

● **IT is 18 colors in all.**

You can mix stainers to have your favorite color.

When you want to weaken a color, mix "Clear" stainer.



※ Even if you apply to white flower, please process from decolorization.



Blue



Red



Yellow



Pink



Green



Purple



Chocolate



Orange



Lavender



Melon



Cobalt



Clear



Bordeaux



SAKURA



Peach

Processing liquid for Preserved flowers.



# Mysty

Palace Chemical CO., LTD  
<http://www.the-flower.biz>